



Product Datasheet

CD303 Antibody (orb1742787)


Description	CD303 is an approximately 38 kDa type II transmembrane glycoprotein with an extracellular C-type lec...
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IHC-Fr
Immunogen	CD303 ectodomain fused with human IgG Fc domain
Target	CD303
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Uniprot ID	Q8WTT0
Entrez	170482
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD303
positive
CD123...



Flow
cytometry
surface
staining
pattern ...