

Product Datasheet

ACE2 (18-740) Protein (orb1742916)

Description	ACE2 (18-740) Protein
Preservatives	Recombinant protein stored in 50mM Tris-HCl, pH 7.5, 90mM glycine.. Recombinant protein stored in 50mM Tris-HCl, pH 7.5, 90mM glycine.
Form/Appearance	liquid
Storage	Stored at -70°C. Stable for 1 year at -70°C.
Tag	Fc
Note	For research use only
Source	rat
MW	110.4
Uniprot ID	Q5EGZ1
Expiration Date	6 months from date of receipt.