

Product Datasheet

2019-nCoV Spike protein S1 (T19R, G142D, del157-158, L452R, T478K, D614G, P681R) (orb1742923)

Description	2019-nCoV Spike protein S1 (T19R, G142D, del157-158, L452R, T478K, D614G, P681R)
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole.
Form/Appearance	liquid
Storage	Stored at -70°C. Stable for 1 year at -70°C.
Tag	His
Note	For research use only
Source	2019-nCoV (virus)
MW	76.2
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.