

Product Datasheet

2019-nCoV spike protein S1 (?Y145, ?H146, E484K, D614G) (orb1742932)

Description	2019-nCoV spike protein S1 (?Y145, ?H146, E484K, D614G)
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole.. Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 150mM imidazole.
Form/Appearance	liquid
Storage	Stored at -70°C. Stable for 1 year at -70°C.
Tag	His
Note	For research use only
Source	2019-nCoV (virus)
MW	76.0
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.