

Product Datasheet

2019-nCoV 3CLpro, Active (orb1742977)

Description	2019-nCoV 3CLpro, Active
Preservatives	Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 200mM imidazole, 1 mM DTT, 10% glycerol.. Recombinant protein stored in 50mM sodium phosphate, pH 7.5, 300mM NaCl, 200mM imidazole, 1 mM DTT, 10% glycerol.
Form/Appearance	liquid
Storage	Stored at -70°C. Stable for 1 year at -70°C.
Tag	His
Note	For research use only
Source	2019-nCoV (virus)
MW	34.8
Expiration Date	6 months from date of receipt.