

Product Datasheet

**Biotinylated Human TIM-3 / HAVCR2 Protein, Fc,
(MALS verified) (orb1742993)**

Description

Biotinylated Human TIM-3

human 293 cells (HEK293). It contains AA Ser...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

TIM-3

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-Fc & C-Avi

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus, followed by an Avi tag. The protein has a calculated MW of 48.0 kDa. The protein migrates as 65-70 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Biotinylated Human TIM-3 Protein, Fc, (orb1742993) is expressed from human 293 cells (HEK293). It contains AA Ser 22 - Arg 200 (Accession # Q8TDQ0-1).

MW

48.0 kDa

Uniprot ID

Q8TDQ0

Expiration Date

6 months from date of receipt.