

## Product Datasheet

**Human MMP-8 Protein, His Tag (active enzyme)**  
**[Out of stock] (orb1742995)**

**Description**

Matrix metalloproteinase-8 (MMP-8)

collagenase and CLG1 is a member of ...

**Reactivity**

Human

**Endotoxins**

1.0 EU per µg

**Conjugation**

Unconjugated

**Target**

MMP-8

**Preservatives**

50 mM Tris, 150 mM NaCl, pH 8.0

**Form/Appearance**

Liquid

**Storage**

-70°C

**Tag**

C-10×His

**Note**

For research use only

**Application notes**

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 29.0 kDa. The protein migrates as 38-55 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Purity**

90%

**Source**

Human MMP-8 Protein, His Tag (orb1742995) is expressed from human 293 cells (HEK293). It contains AA Phe 21 - Gly 262 (Accession # P22894-1).

**MW**

29.0 kDa

**Expiration Date**

6 months from date of receipt.