

Product Datasheet

Human CLEC7A Protein, His Tag (orb1743002)

Description

Human CLEC7A Protein, His Tag is

(HEK293). It contains AA Thr 66 - Me...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CLEC7A

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

N-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 17.5 kDa. The protein migrates as 15-17 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).

Purity

90%

Source

Human CLEC7A Protein, His Tag (orb1743002) is expressed from human 293 cells (HEK293). It contains AA Thr 66 - Met 201 (Accession # NP_072092).

MW

17.5 kDa

Uniprot ID

NP_072092

Expiration Date

6 months from date of receipt.