

Product Datasheet


Cynomolgus IL-12 R beta 1 Protein, Fc Tag (orb1743005)

Description

The human IL-12 R subunit is a member of the

Cynomolgus IL-12 R beta 1 Protein, Fc Tag (83.9 kDa)
1.0 µg of Protein in 200 µL of Sterile PBS, pH 7.4 per vial

	high-a...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-12 R beta 1
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 83.9 kDa. The protein migrates as 110-120 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	95%
Source	Cynomolgus IL-12 R beta 1 Protein, Fc Tag (orb1743005) is expressed from human 293 cells (HEK293). It contains AA Cys 24 - Asp 545 (Accession # G7PWY9-1).
MW	83.9 kDa
Expiration Date	6 months from date of receipt.



Immobilized
Human IL-12B & IL-12A
Hetero...