

Product Datasheet

Human IL-22BP Protein, Fc Tag (MALS verified)
(orb1743040)

Description

This gene encodes a member of the class II

protein spe...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-22BP

Preservatives

PBS, pH 7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;
-70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51 kDa. The protein migrates as 65-85 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source


Human IL-22BP Protein, Fc Tag (orb1743040) is expressed from human 293 cells (HEK293). It contains AA Thr 22 - Pro 231 (Accession # Q969J5-2).

MW

51 kDa

Expiration Date

6 months from date of receipt.



Immobilized
Human IL-
22, His Tag
(Cat. N...