

Product Datasheet

Rhesus macaque CEACAM-6 / CD66c Protein, His Tag (orb1743050)

Description

Rhesus macaque CEACAM-6, His Tag

Reactivity

Monkey

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CEACAM-6

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 33.5 kDa. The protein migrates as 50-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

XP_014979566.2

Purity

90%

Source

Rhesus macaque CEACAM-6, His Tag (orb1743050) is expressed from human 293 cells (HEK293). It contains AA Gln 35 - Gly 320 (Accession # XP_014979566.2).

MW

33.5 kDa

NCBI
[XP_014979566.2](#)
Expiration Date

6 months from date of receipt.