

Product Datasheet

Human MRGPRX2 full length protein-synthetic nanodisc (orb1743166)

Description	Human MRGPRX2 full length protein-synthetic nanodisc
Reactivity	Human
Target	MRGPRX2
Form/Appearance	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Source	Mammalian
MW	The human full length MRGPRX2 protein has a MW of 37.1 kDa
Uniprot ID	Q96LB1
Expiration Date	6 months from date of receipt.