

Product Datasheet

Human CCR7 full length protein-synthetic nanodisc (orb1743211)

Description	Human CCR7 full length protein-synthetic nanodisc
Reactivity	Human
Target	CCR7
Form/Appearance	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Source	Mammalian
MW	The human full length CCR7 protein has a MW of 42.9 kDa
Uniprot ID	P32248
Expiration Date	6 months from date of receipt.