

Product Datasheet


Human CCR8 full length protein-synthetic nanodisc (orb1743234)

Description


Human CCR8 full length protein-synthetic

ELISA assay to evaluate CCR8 Nanodisc
Flag tagged CCR8 Nanodisc per well


Reactivity	Human
Target	CCR8
Form/Appearance	Lyophilized from nanodisc solubilization buffer (20 mM Tris-HCl, 150 mM NaCl, pH 8.0). Normally 5% - 8% trehalose is added as protectants before lyophilization
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing)
Note	For research use only
Source	HEK293
MW	The human full length CCR8 Protein has a MW of 40.7 kDa
Uniprot ID	P51685
Expiration Date	6 months from date of receipt.



Elisa plates were pre-coated with Flag T...



WB analysis of Human CCR8-Nanodisc with ...



Human CCR8-Nanodisc can bind Anti-CCR8 a...