

Product Datasheet

ERBB2 Antibody (orb1743595)

Description	ERBB2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	A synthesized peptide derived from Human ERBB2
Target	ERBB2
Preservatives	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Mouse IgG
Clonality	Monoclonal
Uniprot ID	P04626
Expiration Date	12 months from date of receipt.