



Product Datasheet

Anti-F4/80/Adgre1 Antibody (orb1743751)


Description	Anti-F4/80/Adgre1 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	E.coli-derived mouse F4/80/Adgre1 recombinant protein (Position: G27-S230).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml
Clonality	Polyclonal
Uniprot ID	Q61549
Dilution Range	Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Mouse, Rat Immunofluorescence, 5 µg/ml, Mouse, Rat
Expiration Date	12 months from date of receipt.



IF analysis of F4/80/Adgre1 using anti-F...



IF analysis of F4/80/Adgre1 using anti-F...



IF analysis of F4/80/Adgre1 using anti-F...