

Product Datasheet

Anti-NUP37 Antibody (orb1743787)

Description

Anti-NUP37 Antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested Applications

ELISA, FC, ICC, IF, IHC, WB

Immunogen

E.coli-derived human NUP37 recombinant protein (Position: R7-D231).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml

Clonality

Polyclonal

MW

37 kDa

Uniprot ID


[Q8NFH4](#)

Dilution Range

Western blot, 0.25-0.5 µg/ml, Human
Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human, Rat
Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human
Flow Cytometry (Fixed), 1-3 µg/1x10⁶ cells, Human Direct
ELISA, 0.1-0.5 µg/ml, Human

Expiration Date

12 months from date of receipt.



Flow Cytometry analysis of MCF-7 cells

