


Product Datasheet

Human MSR ELISA Kit (orb1743979)

Description	Human MSR ELISA Kit (96 Tests). Quantitate Human MSR in cell culture supernatants, serum and plasma ...
Reactivity	Human
Range	93.7 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K77-L451
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	P21757
Sensitivity	45 pg/ml
Expiration Date	Please enquire.



Human MSR
PicoKine
ELISA Kit
standard
CU...