


Product Datasheet

Human Napsin A ELISA Kit (orb1743984)

Description	Human Napsin A ELISA Kit (96 Tests). Quantitate Human Napsin A in cell culture supernatants, serum a...
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: T25-G420
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin)
Uniprot ID	O96009
Sensitivity	30 pg/ml
Expiration Date	Please enquire.



Human
Napsin A
PicoKine
ELISA Kit
standa...