

Product Datasheet

CGRP Antibody (orb175801)

Description	CGRP Antibody
Species/Host	Gallus
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	KDNFVPTNVGSEAF-NH2
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Chicken IgY
Uniprot ID	075575
Entrez	27297
Expiration Date	12 months from date of receipt.