

Product Datasheet

GFAP Antibody (orb175805)

Description	GFAP Antibody
Species/Host	Gallus
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	Purified bovine spinal cord GFAP
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Chicken IgY
Uniprot ID	P14136
Entrez	2670
Dilution Range	ICC 1:2000; IHC 1:2000; WB 1:10,000
Expiration Date	12 months from date of receipt.