

Product Datasheet

Neurofilament NF-M Antibody (orb175809)

Description	Neurofilament NF-M Antibody
Species/Host	Gallus
Reactivity	Gallus, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	C-terminal extension of rat NF-M, the so-called KE segment, was expressed in bacteria and purified.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Chicken IgY
Uniprot ID	P07197
Entrez	4741
Dilution Range	ICC 1:2000; IHC 1:2000; WB 1:10,000
Expiration Date	12 months from date of receipt.