



Product Datasheet

Vimentin Antibody (orb175810)


Description	Vimentin Antibody
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant vimentin purified from E. Coli
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Chicken IgY
Uniprot ID	M14144
Entrez	7431
Expiration Date	12 months from date of receipt.



WB of crude extract of human embryonic k...



hN2 Primary Neurons stained with our chi...



Human HeLa cells stained with Vimentin a...