

## Product Datasheet

### proDynorphin Antibody (orb175853)

|                            |                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>         | proDynorphin Antibody                                                                    |
| <b>Species/Host</b>        | Guinea pig                                                                               |
| <b>Reactivity</b>          | Mouse, Rat                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                             |
| <b>Tested Applications</b> | IHC                                                                                      |
| <b>Form/Appearance</b>     | Affinity Purified                                                                        |
| <b>Storage</b>             | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                    |
| <b>Isotype</b>             | Guinea Pig IgG                                                                           |
| <b>Uniprot ID</b>          | <b>P01213</b>                                                                            |
| <b>Entrez</b>              | <b>5173</b>                                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                          |