

Product Datasheet

ppNPY Antibody (orb175861)

Description	ppNPY Antibody
Species/Host	Guinea pig
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Form/Appearance	Protein G Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Guinea Pig IgG
Uniprot ID	P01303
Entrez	4852
Expiration Date	12 months from date of receipt.