

Product Datasheet

Synaptophysin Antibody (orb175877)

Description	Synaptophysin Antibody
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synaptophysin c-terminal peptide
Form/Appearance	Supernatant
Storage	Store at 4°C short term. Store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	P08247
Entrez	6855
Expiration Date	12 months from date of receipt.