

Product Datasheet

Calbindin Antibody (orb175890)

Description	Calbindin Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Prokaryotic recombinant protein corresponding to the majority of the calbindin molecule.
Form/Appearance	Supernatant
Storage	Store at 4°C short term. Store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	P05937
Entrez	793
Expiration Date	12 months from date of receipt.