

Product Datasheet

ChAT Antibody (orb175893)

Description	ChAT Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Prokaryotic recombinant protein corresponding to the C-terminal region of the human ChAT molecule.
Form/Appearance	Supernatant
Storage	Store at 4°C short term. Store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	P28329
Entrez	1103
Dilution Range	ICC 1:100-1:200; IHC 1:100-1:200; WB 1:250-1:500
Expiration Date	12 months from date of receipt.