

Product Datasheet

Peripherin Antibody (orb175916)

Description	Peripherin Antibody
Species/Host	Mouse
Reactivity	Feline, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant rat peripherin purified from E. coli. .Clone: 7C5.
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG
Uniprot ID	P41219
Entrez	5630
Dilution Range	ICC 1:200; IHC1:200; WB 1:1000
Expiration Date	12 months from date of receipt.