

## Product Datasheet

### VR1 peptides (orb175989)

|                        |                                                                                          |
|------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>     | VR1 peptides                                                                             |
| <b>Form/Appearance</b> | Liquid                                                                                   |
| <b>Storage</b>         | Store at 4°C short term. Store at -20°C for maximum stability. Avoid freeze-thaw cycles. |
| <b>Note</b>            | For research use only                                                                    |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                          |