

Product Datasheet

Bradykinin Receptor B2 Antibody (orb176136)

Description	Bradykinin Receptor B2 Antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	CMGESVQMENSMGTLR
Form/Appearance	Whole serum
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P30411
Entrez	624
Expiration Date	12 months from date of receipt.