

Product Datasheet

Humanin Antibody (orb176170)

Description	Humanin Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	LLLLTSEIDLPAKRRR (a.a. 9-24 of Humanin peptide)
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	Q8IVG9
Entrez	4550
Expiration Date	12 months from date of receipt.