

Product Datasheet

CLAC-P (NC2-1 region) Antibody (orb176173)

Description	CLAC-P (NC2-1 region) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	NHGFLSADQQLIK (a.a. 171-183); NC2-1 region of human and mouse CLAC-P
Form/Appearance	Affinity Purified
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	Q9BXS0
Entrez	84570
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.