

Product Datasheet

Peripherin Antibody (orb176252)

Description	Peripherin Antibody
Species/Host	Rabbit
Reactivity	Feline, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant full length peripherin purified from E. coli.
Form/Appearance	Whole serum
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P41219
Entrez	5630
Dilution Range	ICC 1:100-1:200; IHC 1:100-1:200; WB 1:250-1:500
Expiration Date	12 months from date of receipt.