

Product Datasheet

MARCKS (Myristoylated Alanine Rich C Kinase Substrate) Antibody (orb176254)

Description	MARCKS (Myristoylated Alanine Rich C Kinase Substrate) Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Purified recombinant MARCKS protein.
Form/Appearance	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P29966
Entrez	4082
Dilution Range	ICC 1:1000; IHC 1:1000; WB 1:10,000
Expiration Date	12 months from date of receipt.