

Product Datasheet

CGRP Antibody (orb176264)

Description	CGRP Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC
Immunogen	Synthetic rat alpha-CGRP coupled to bovine thyroglobulin (BTg) with glutaraldehyde.
Form/Appearance	Whole serum
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P01258
Entrez	796
Expiration Date	12 months from date of receipt.