

Product Datasheet

EGR2 Antibody (orb176298)

Description	EGR2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide made to a portion of human EGR2 protein sequence within residues 200-300. [Swiss-Prot# P11161].
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P11161
Entrez	1959
Expiration Date	12 months from date of receipt.