



## Product Datasheet

### SSEA3 Antibody (orb176301)

|                            |                                                       |
|----------------------------|-------------------------------------------------------|
| <b>Description</b>         | SSEA3 Antibody                                        |
| <b>Species/Host</b>        | Rabbit                                                |
| <b>Reactivity</b>          | Human, Mouse                                          |
| <b>Conjugation</b>         | Unconjugated                                          |
| <b>Tested Applications</b> | FC, ICC, IP, WB                                       |
| <b>Immunogen</b>           | 4-8 cell stage mouse embryos.                         |
| <b>Form/Appearance</b>     | liquid                                                |
| <b>Storage</b>             | Aliquot and store at -20°C. Avoid freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                 |
| <b>Isotype</b>             | Rabbit IgG                                            |
| <b>Expiration Date</b>     | 12 months from date of receipt.                       |



Immunofluorescence analysis of HeLa cell...



Flow cytometric analysis of NTERA-2 cel...