

Product Datasheet

FE65 Antibody (orb176302)

Description	FE65 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide corresponding to residues surrounding Thr668 of human APP695
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	000213
Entrez	322
Expiration Date	12 months from date of receipt.