

Product Datasheet

TRPM2 Antibody (orb176304)

Description	TRPM2 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide made to a C-terminal portion of the rat protein (within residues 1430-1508). [Swiss-Prot# Q5G856].
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	094759
Entrez	7226
Expiration Date	12 months from date of receipt.