

Product Datasheet

CD36 Antibody (orb176305)

Description	CD36 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic peptide mapping to a region of human CD36 between residues 100-200.
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P16671
Entrez	948
Expiration Date	12 months from date of receipt.