

Product Datasheet

NMDA Receptor 2A Antibody (orb176307)

Description	NMDA Receptor 2A Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	C-terminal fusion protein representing amino acids 1253 - 1391 of rat NMDAR2A.
Form/Appearance	Whole serum
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	Q12879
Entrez	2903
Expiration Date	12 months from date of receipt.