

Product Datasheet

LC3 Antibody (orb176324)

Description	LC3 Antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide made to an internal portion of the human LC3 protein sequence (between residues 50-125).
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	Q9H492
Entrez	84557
Expiration Date	12 months from date of receipt.