

Product Datasheet

GLUT1 Antibody (orb176339)

Description	GLUT1 Antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Primate, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide to an internal region (within residues 300-400) of the human LPIN1 protein. [Swiss-Prot# Q14693].
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P11166
Entrez	6513
Expiration Date	12 months from date of receipt.