

Product Datasheet

Olig2 Antibody (orb176341)

Description	Olig2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	Recombinant mouse Olig2.
Form/Appearance	Protein A purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	Q13516
Entrez	10215
Expiration Date	12 months from date of receipt.