

Product Datasheet

Aggrecan (ACAN) Neoepitope Antibody (orb176342)

Description	Aggrecan (ACAN) Neoepitope Antibody
Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic peptide corresponding to residues C G G N(a.a. 387) I T E G E(a.a. 392) of human aggrecan.
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Uniprot ID	P16112
Entrez	176
Expiration Date	12 months from date of receipt.