

Product Datasheet

Tyrosine Hydroxylase Antibody (orb176353)

Description	Tyrosine Hydroxylase Antibody
Species/Host	Sheep
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Native rat tyrosine hydroxylase purified from pheochromocytoma.
Form/Appearance	Affinity Purified
Storage	Aliquot and store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Isotype	Sheep IgG
Uniprot ID	P07101
Entrez	7054
Dilution Range	ICC 1:100-1:200; IHC 1:100-1:200; WB 1:250-1:500
Expiration Date	12 months from date of receipt.