



Product Datasheet

ADAM19 antibody (orb176701)

Description	Goat polyclonal antibody to ADAM19.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ADAM19
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 125kDa band observed in Human Spleen lysates (calculated MW of 101kDa according to NP_150377.1). The observed molecular weight corresponds to the glycosylated form. Recommended concentration: 0.1-0.3µg/ml.. IHC: Paraffin embedded Human Spleen. Recommended concentration: 5µg/ml.
Protein Sequence	KIQQSSEARPLESN
Clonality	Polyclonal
MW	101
Entrez	8728
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Spleen lysate u...



Immunohistochemical staining of Human Sp...