



Product Datasheet

CLIC1 antibody (orb176709)

Description	Goat polyclonal antibody to CLIC1.
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	CLIC1 / NCC27
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 30kDa band observed in lysates of cell line U937 (calculated MW of 26.9kDa according to NP_001279.2). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Colon. Recommended concentration: 5µg/ml.
Protein Sequence	TEVHTDTNKIEE
Clonality	Polyclonal
MW	26.9
Entrez	1192
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of U937 lysate using CLIC1 antibody.



Immunohistochemical staining of Human Colon tissue using CLIC1 antibody.