



Product Datasheet

cytokeratin 19 antibody (orb176715)

Description	Goat polyclonal antibody to cytokeratin 19.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	cytokeratin 19 (aa285-298)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 43kDa band observed in Human cancer lysates of Breast, Pancreas and Colon (calculated MW of 44.1kDa according to NP_002267.2). Recommended concentration: 0.03-0.1µg/ml. Primary incubation was 1 hour.. IHC: Paraffin embedded Human Liver. Recommended concentration: 2.5µg/ml.
Protein Sequence	HTEQLQMSRSEVTD
Clonality	Polyclonal
MW	44.1
Entrez	3880
Dilution Range	ELISA: 1:128000, WB: 0.03-0.1 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis
of Human Colon
can...



Immunohistochemical
staining of Human
Li...